Semi-automatic Coagulometer

Coagulometer 2

Versatile coagulometer controlled by microprocessor with turbodensidometric detection of the beginning of the Clot formation (fibrin polymerization).

The clot detection is carried out in positive O.D. variations or even in negative ones in high level fibrinogenemia.

It is very easy to use and can store up to 10 methods in memory.

Performance of all commonly used tests like Quick test (PT), APTT, Fibrinogen, TT, ATIII, factor assays, etc., with a high level of precision and repeatability.

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TWO CHANNELS COAGULOMETER

The detection lamp intensity is automatically adjusted on each channel after addition of the initial reagent on the measuring cell. Consequently the baseline is adjusted as well. Automatic Clot detection is carried out in a homogeneous mixed medium by means of a system that does not interfere with the normal clot formation and avoids possible sedimentation phenomena of activating reagents during the test, with the turbodetection concentrated in the central axis of the cuvette.

All results are displayed and automatically printed out on continuous paper, together with the ID number, the used method and the alarms for each channel. Results can be given one by one, or as a mean value between two consecutive readings.

TECHNICAL SPECIFICATIONS

Optical group: The optical group of each channel is composed by a light source with automatic adjustment of its intensity, a diffusion lens, a diaphragm to control the light path, a reaction cell and a high sensitivity revealing element which detects all O.D. variations in the centre of the cuvette.

Thermostatic Block: Thermostated device at $37 \pm 0,1^\circ\text{C}$ that includes the two optical systems and the two magnetic stirrer motors. It allows to thermostat up to 21 cuvettes and two vials of reagent.

Numeric keyboard: Easy to handle, to program and to perform the tests.

Alphanumeric display: LCD of 4 x 40 characters in which all results, alarm indications, temperature control, control times and all data of the test to be done, appear.

Printer: Thermal using 56 mm thermal paper.

Storage capacity: Up to 10 methods.

Program capacity: Incubation time, delay time, security time and calibration curves up to nine points for the methods that require it.

Calculations: It automatically calculates results related to the stored calibration curves with the corresponding units, activity and International Normalized Ratio (INR).

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