



AGGLUTINOSCOPE

MOD. 705

COD. 29980001



This equipment has been designed, built and tested following the European Directive 2014/30/EU, 2014/35/EU, 2006/42/CE, 2011/65/EU, the normative EN 61010-1:2010, EN 61326-1:2013, EN ISO 12100:2010, with the following National Laws D.lgs 81/08 del 09/04/2008, D.lgs 14 marzo 2014 n. 49

DESCRIPTION

The Agglutinoscope model 705 consists of a tilting body made in powdered painted steel acids resistant with an acrylic resin plate to provide uniform diffusion of light.

With the regulator is possible to settle luminosity and temperature.

Thermometer adjusted to 37°C.

Main switch bright.

FEATURES

TECHNICAL SPECIFICATIONS

- ♦ Work surface dimensions WxD: 250 x 85 mm.
- ♦ Overall dimensions WxDxH: 360 x 180 x 150 mm.
- ♦ Weight: 3 Kg

ELECTRICAL SPECIFICATIONS

- ♦ Power supply: 230 V - 50 Hz.
- ♦ Electrical input: 40 W (n°2 lamps 2x18 W)
- ♦ Overload fuses: 2 x 0.5 A F (5x20) mm.
- ♦ Connection outlet: 10 A

