

**MAXI MIXER VORTEX**

**MOD. 714/R**

**WITH SPEED REGULATOR**

**COD. 29970030**



CE

<b>DESCRIPTION</b>
--------------------

The Maxi Mixer Vortex mod.714/R shakes test tubes of different shapes and dimensions.

Equipped with speed regulator from 50 to 900 rpm.

Outer casing in powdered painted steel acids resistant.

Neoprene shaking cup Ø 25 mm. Rubber suction pads.

Eccentricity vibration 3.7 mm.

Main switch bright.

### **TECHNICAL SPECIFICATIONS**

- ♦ Overall dimensions WxDxH: 130 mm x 200 x 150 mm.
- ♦ Weight: 4 Kg
- ♦ Shaking cup speed: from 50 to 900 rpm

### **ELECTRICAL SPECIFICATIONS**

- ♦ Power supply: 230 V - 50 Hz.
- ♦ Electrical input: 70 W
- ♦ Overload fuses: 2 x 0.5 AF (5 x 20) mm.
- ♦ Connection outlet: 10 A

Neoprene shaking cup Ø25 mm



Speed regulator

Overload fuses

Main switch bright